

# HEATING MANTLES

- » for operating temperatures up to 450°C
- » for glass flasks, round bottom up to 20 l.
- » flexible glass yarn heating element
- » plastic coated metal housing
- » thermally insulated and grounded
- » predrilled threaded hole on rear allows connection of the clamp KM-SK (accessory)

## BOECO **LabHEAT**<sup>®</sup> KM-M SERIES

1,5 m power supply (earthed) cable with heating-zone switch



KM-M

Volume ml/l	Flask ø mm	Power W	Heating zones	Order No. 220 V AC	Order No. 110 V AC
50	51	55	1	BOE 91000050	BOE 91500050
100	64	100	1	BOE 91000100	BOE 91500100
250	85	150	2	BOE 91000250	BOE 91500250
500	105	200	2	BOE 91000500	BOE 91500500
1	131	300	2	BOE 91001000	BOE 91501000
2	166	500	2	BOE 91002000	BOE 91502000
3	185	600	2	BOE 91003000	BOE 91503000
4	207	750	2	BOE 91004000	BOE 91504000
5	223	860	2	BOE 91005000	BOE 91505000
6	236	1000	2	BOE 91006000	BOE 91506000
10	279	1400	2	BOE 91010000	BOE 91510000
20	345	2000	2	BOE 91020000	BOE 91520000

## BOECO **LabHEAT**<sup>®</sup> KM-ME SERIES

designed identical to series KM-M, but in addition equipped with a built-in heating-zone switch and a power controller which allows a continuous adjustment of the heater power.



KM-ME

Volume ml/l	Flask ø mm	Power W	Heating zones	Order No. 220 V AC	Order No. 110 V AC
50	51	55	1	BOE 92000050	BOE 92500050
100	64	100	1	BOE 92000100	BOE 92500100
250	85	150	2	BOE 92000250	BOE 92500250
500	105	200	2	BOE 92000500	BOE 92500500
1	131	300	2	BOE 92001000	BOE 92501000
2	166	500	2	BOE 92002000	BOE 92502000
3	185	600	2	BOE 92003000	BOE 92503000
4	207	750	2	BOE 92004000	BOE 92504000
5	223	860	2	BOE 92005000	BOE 92505000
6	236	1000	2	BOE 92006000	BOE 92506000
10	279	1400	2	BOE 92010000	BOE 92510000
20	345	2000	2	BOE 92020000	BOE 92520000

## BOECO **LabHEAT**<sup>®</sup> KM-MER SERIES

designed identical to series KM-ME, but in addition equipped with a built-in separate switched magnetic stirrer with an infinitely variable speed range - up to 1600 rpm.



KM-MER

Volume ml/l	Flask ø mm	Power W	Heating zones	Order No. 220 V AC
100	64	100	1	BOE 94000100
250	85	150	2	BOE 94000250
500	105	200	2	BOE 94000500
1	131	300	2	BOE 94001000



Code	Description
BOE 970000001	Support clamp KM-SK for the attachment of the metal-cased heating mantles to support rod or within support wall.

## BOECO

The family-owned company Boeckel + Co. (GmbH + Co.) KG was founded in 1929 and now exists in the fourth generation.

To our customers on all 5 continents we are pleased to introduce our fifth general BOECO catalog with our extended and updated complete line of laboratory equipment and consumables.

We offer high quality at competitive prices by a convenient and cost-effective one-stop supply.

## CONTACT

**Boeckel + Co. (GmbH + Co.) KG**  
**Roedingsmarkt 33**  
**20459 Hamburg**  
**Germany**

**Phone:** + 49 40 325627-0  
**Fax:** + 49 40 325627-41  
**E-mail:** [info@boeco.com](mailto:info@boeco.com)  
**Internet:** [www.boeco.com](http://www.boeco.com)

**V.A.T. ID-No: DE 118903538**

## BANK

### EURO €

Deutsche Bank AG, Hamburg,  
SWIFT: DEUTDEHHXXX  
IBAN- No. DE 65 2007 0000 0010 5601 00

### USD \$

Deutsche Bank AG, Hamburg,  
SWIFT: DEUTDEHHXXX  
IBAN- No. DE 27 2007 0000 0010 5601 05